Nicolai Welfler, Product manager

Vitalthings

- et norsk helseteknologiselskap

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¹Ståle Tøften, Ståle Pallesen, Maria Hrozanova, Frode Moen, Janne Grønli: "Validation of sleep stage classification using noncontact radar technology and machine learning (Somnofy®)", available online 6 March 2020, epub ahead of print, Sleep Medicine. DOI: 10.1016/j.sleep.2020.02.022

Original Article

Validation of sleep stage classification using non-contact radar technology and machine learning (Somnofy®)

Ståle Tøften ^{a,*}, Ståle Pallesen ^{b,c}, Maria Hrozanova ^d, Frode Moen ^e, Janne Grønli ^f

^a Department of Data Science, VitalThings AS, Tønsberg, Norway

^b Department of Psychosocial Science, University of Bergen, Bergen, Norway

^c Norwegian Competence Center for Sleep Disorders, Haukeland University Hospital, Bergen, Norway

^d Department of Neuromedicine and Movement Science, Faculty of Medicine and Health Sciences, Centre of Elite Sport Research, NTNU, Trondheim, Norway

^e Department of Education and Lifelong Learning, Centre of Elite Sport Research, NTNU, Trondheim, Norway

^f Department of Biological and Medical Psychology, University of Bergen, Bergen, Norway

ARTICLE INFO

Article history:

Received 8 November 2019

Received in revised form

11 February 2020

Accepted 23 February 2020

Available online 6 March 2020

Keywords:

Polysomnography

Radar

Sleep stages

Somnofy

Validation

ABSTRACT

Objective: To validate automatic sleep stage classification using deep neural networks on sleep assessed by radar technology in the commercially available sleep assistant Somnofy® against polysomnography (PSG).

Methods: Seventy-one nights of overnight sleep in healthy individuals were assessed by both PSG and Somnofy at two different institutions. The Somnofy unit was placed in two different locations per room (nightstand and wall). The sleep algorithm was validated against PSG using a 25-fold cross validation technique, and performance was compared to the inter-rater reliability between the PSG sleep scored by two independent sleep specialists.

Results: Epoch-by-epoch analyses showed a sensitivity (accuracy to detect sleep) and specificity (accuracy to detect wake) for Somnofy of 0.97 and 0.72 respectively, compared to 0.99 and 0.85 for the PSG scorers. The sleep stage differentiation for Somnofy was 0.75 for N1/N2, 0.74 for N3 and 0.78 for R, whilst PSG scorers ranged between 0.83 and 0.96. The intraclass correlation coefficient revealed excellent and good reliability for total sleep time and sleep efficiency, while sleep onset and R latency had poor agreement. Somnofy underestimated total wake time by 5 min and N1/N2 by 3 min. N3 was overestimated by 4 min and R by 3 min. Results were independent of institution and sensor location.

Conclusion: Somnofy showed a high accuracy staging sleep in healthy individuals and has potential to assess sleep quality and quantity in a sample of healthy, mostly young adults. More research is needed to examine performance in children, older individuals and those with sleep disorders.

© 2020 Elsevier B.V. All rights reserved.

1. Introduction

The gold standard for objective sleep assessment, polysomnography (PSG), is typically performed in a sleep laboratory and data needs to be manually scored by a sleep technician, making it obtrusive, costly and less suitable for longitudinal studies. Hence, there is a need for validated low-cost equipment for the assessment of sleep, which is user-friendly, accurate and non-intrusive. From a clinical and research perspective, the capacity to obtain longitudinal sleep–wake data may individualize treatment decision and health optimization and improve disease phenotyping.

Radar technology has a great potential for home sleep assessment as it is completely non-intrusive. Technology based on impulse radio ultra-wideband (IR-UWB) radar sensor technique has been shown to reliably monitor vital signs in real-time as respiratory and cardiac events, as well as limb movements [1–4]. Respiration and movement have been shown to correlate with REM and non-REM sleep [5,6]. Recently, an IR-UWB radar was developed for the purposes of sleep assessment [7]. In the pilot validation study against PSG, the IR-UWB radar quantified sleep and wakefulness by an algorithm integrating movements from all body parts. Results of the study revealed a small overall discrepancy between PSG estimates for total sleep time, and the mean sensitivity (radar = sleep when PSG = sleep) and specificity (radar = wake when PSG = wake) were higher or comparable to that reported for actigraphic studies [8]. Despite the promising results of the IR-UWB

* Corresponding author. Department of Data Science, VitalThings AS, Jarlsveien 48, Tønsberg, Vestfold 3124, Norway.
E-mail address: st@vitalthings.com (S. Tøften).

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